

For safety reasons the power cord should be kept under the bed or away from foot traffic.

CPM Machine should be kept off the floor for infection control purposes.

Follow physicians order on frequency and duration.

EXTENSION



can be changed by Pressing and Holding the Extension button and then pressing the up or down button to increase or decrease to the desired degree.

FLEXION



can be changed by Pressing and Holding the Flexion button and then pressing the up or down button to increase or decrease to the desired degree.

CARRIAGE SPEED



can be changed by Pressing and Holding the Carriage Speed button and then pressing the up or down button to increase or decrease to the desired degree.

START/STOP



This button functions as BOTH the START and STOP switch.

Pressing this button once will STOP the unit and pressing this button a second time will REVERSE the direction of the carriage.

The machine pads are single patient use and must be discarded at the completion of your therapy.

The use of the CPM machine is commonly used from 10 to 21 days depending on your insurance carrier.

I understand that upon the completion of authorized length of time, If I so desire to continue with the therapy I must contact my insurance carrier for an extension of time.